



ACTIVELY DELIVERING HIGH PERFORMANCE WOUND CARE SOLUTIONS

Aquafiber Extra is 68% stronger and has superior fluid handling capabilities compared to Aquacel Extra.*



activheal.com
AQUAFIBER EXTRA

*Based on in vitro testing, AMS data on file, LD030-16, LD024-16, LD037-16, LD093-15

HIGHER FLUID RETENTION

The unique design & fibre engineering of ActivHeal® Aquafiber Extra allows the dressing to absorb more exudate than Aquacel Extra, transforming into a cohesive gel, locking fluid away from the wound.

AQUACEL
EXTRA
19.8
g/100cm²

AQUAFIBER
EXTRA
25.9
g/100cm²

AQUAFIBER
EXTRA
7.6
N/cm

AQUACEL
EXTRA
4.5
N/cm

EVEN STRONGER WHEN WET

The quilting pattern of ActivHeal® Aquafiber Extra result in a 68% stronger dressing than Aquacel Extra.

HIGHER ABSORBENCY

ActivHeal® Aquafiber Extra has higher absorbency compared to other leading competitors

AQUAFIBER
EXTRA
36.4
g/100cm²


KERRACEL
30.6
g/100cm²

AQUACEL
EXTRA
24.7
g/100cm²


DURAFIBER
24.6
g/100cm²

AMS data on file LD030-16, LD024-16, LD037-16, LD09-15

NHSSC (Inc VAT)



DRUG TARIFF



10cm x 10cm

Advanced Medical Solutions Ltd

Premier Park, 33 Road One, Winsford Industrial Estate, Winsford, Cheshire CW7 3RT
T 08444 125 755 E activheal@admedsol.com W www.activheal.com

